

**Antibody to Dynein Light Chain 1/3 (10F6-D1-D9)**  
MOUSE MONOCLONAL

**Catalog Number:** AB-V73  
**Quantity:** 100 micrograms  
**Format:** PBS (0.14 M Sodium Chloride; 0.003 M Potassium Chloride; 0.002 M Potassium Phosphate; 0.01 M Sodium Phosphate; pH 7.4), no preservative.  
**Host:** Mouse  
**Isotype:** IgG1  
**Clone:** 10F6-D1-D9  
**Immunogen:** Dynein light chain

**Background:** Eukaryotic cells depend on actin- and microtubule-based motor proteins to drive intracellular transport. Dyneins are large multi-subunit protein complexes that move cargo toward the minus ends of microtubules. Cytoplasmic dynein plays critical roles in transporting membrane-bound organelles such as endosomes, lysosomes, and mitochondria, as well as in positioning the centrosome, nucleus, and Golgi apparatus. It is also involved in virus transport to the nucleus, retrograde axonal transport, and the movement of microtubules and neurofilaments. Three families of dynein light chains (DYNL) have been described, with DYNLT family members DYNLT1 and DYNLT3 proposed to mediate the attachment of dynein to specific cargoes.

**Specificity & Preparation:** This IgG1 mouse antibody is generated against Dynein Light Chain and recognizes rat and human homologs of both DYNLT1 and DYNLT3.

**Usage:** Applications include immunoprecipitation and immunofluorescence. Working dilutions must be determined by end user.

**Storage:** Store antibody at -20°C for one year. Avoid repeated freezing and thawing. Gently spin down material 5-10 seconds in a microfuge before use.

**Selected References:**

1. Lo KW, Kogoy JM, Rasoul BA, King SM, Pfister KK (2007) Interaction of the DYNLT (TCTEX1/RP3) light chains and the intermediate chains reveals novel intersubunit regulation during assembly of the dynein complex. *J Biol Chem* 282(51):36871-36878. doi: 10.1074/jbc.M705991200 PMID: 17965411

Scan to view  
all product  
references.

To view protocol(s) for this and other products please visit: [www.ATSBio.com/library/protocols](http://www.ATSBio.com/library/protocols)