

**Antibody to Complement C3 (6C9)**
MOUSE MONOCLONAL

Catalog Number: AB-V85
Quantity: 100 micrograms
Format: PBS (0.14 M Sodium Chloride; 0.003 M Potassium Chloride; 0.002 M Potassium Phosphate; 0.01 M Sodium Phosphate; pH 7.4), no preservative.
Host: Mouse
Isotype: IgG1
Clone: 6C9
Immunogen: Human C3

Background: The complement system is a network of small proteins in the blood, primarily synthesized by the liver, that circulate as inactive precursors until activated by specific triggers. Proteolytic activation initiates a cascade that amplifies the immune response and culminates in formation of the membrane attack complex, which can lyse target cells. Complement component 3 (C3) is central to this process. Cleavage of C3 by C3 convertase generates C3b and iC3b, with C3b driving inflammation and covalently binding to microbial surfaces to promote opsonization.

Specificity & Preparation: This mouse IgG1 monoclonal antibody is generated against human complement component 3 (C3) and recognizes mouse and human homologs of C3 and breakdown products C3b, iC3b and C3d. 6C9 recognizes human & mouse C3b and iC3b.

Usage: Applications include western blot, immunoprecipitation, flow cytometry, and ELISA. Working dilutions must be determined by end user.

Storage: Store antibody at -20°C for one year. Avoid repeated freezing and thawing. Gently spin down material 5-10 seconds in a microfuge before use.

**Selected References:**

1. Li Y, Williams ME, Cousar JB, Pawluczkoysz AW, Lindorfer MA, Taylor RP (2007) Rituximab-CD20 complexes are shaved from Z138 mantle cell lymphoma cells in intravenous and subcutaneous SCID mouse models. *J Immunol* 179(6):4263-4271. doi: 10.4049/jimmunol.179.6.4263 PMID: 17785867

Scan to view
all product
references.

To view protocol(s) for this and other products please visit: www.ATSBio.com/library/protocols