

**FITC-labeled Goat Anti-Rat IgG (H+L)**
FLUORESCENT CONJUGATE

Catalog Number: FL-10
Quantity: 500 micrograms
Format: 50% PBS (0.14 M Sodium Chloride; 0.003 M Potassium Chloride; 0.002 M Potassium Phosphate; 0.01 M Sodium Phosphate; pH 7.4), 50% glycerol; no preservative.
Host: Goat

Background: FITC-labeled goat anti-rat IgG can be used to verify specific binding of rat IgG to its receptor. By first incubating cells with the primary rat antibody, and then binding the FITC-labeled goat anti-rat IgG to the primary antibody, a fluorescent marker is formed that can demonstrate expression of a receptor or affinity of an antibody for its receptor. FITC is excited by 488 nm wavelength light, and emits at 525 nm.

Specificity & Preparation: This fluorescent conjugate was prepared using goat anti-rat IgG (H+L) and the fluorescent compound, fluorescein isothiocyanate (FITC). The antibody binds to rat IgG, and is affinity-purified to decrease background and non-specific binding. This antibody exhibits maximal binding to rat IgG antibodies, and minimal cross-reactivity with other molecules. This product is routinely tested by flow cytometry.

Usage: Applications include flow cytometry (ATS in-house; 2 μ g/ 10^6 cells per 200 μ l).

Storage: Gently spin down material 5-10 seconds in a microfuge before use. The material can be handled safely using normal laboratory precautions. Store the antibody at -20°C for up to one year.

To view protocol(s) for this and other products please visit: www.ATSBio.com/library/protocols