

BIgG-SAP Mouse
NON-TARGETED SAPORIN CONTROL MOLECULE
a tool for use as control for bonded mouse antibody-based lesioning agent

Catalog Number: IT-74
Quantity: 25 micrograms, 100 micrograms, 250 micrograms
Format: PBS (0.14 M Sodium Chloride; 0.003 M Potassium Chloride; 0.002 M Potassium Phosphate; 0.01 M Sodium Phosphate; pH 7.4), no preservative. Sterile-filtered.
Host: Mouse

Background:

Controls are a vital part of the scientific procedure; without them it is difficult to isolate the specific effects from the non-specific or artifactual. This control molecule is the same molecular weight, consists of similar, comparable materials and is synthesized with the same protocols as the targeted conjugates. The difference is the cell-specific targeting agents are replaced with "blanks," antibodies or peptides that have no specificity, and no ability to target cells. In short, they are the perfect control molecules for behavioral experiments with Advanced Targeting Systems' targeted conjugates.

BIgG-SAP mouse is a biotinylated normal mouse IgG mixed with Streptavidin-ZAP in the identical manner as other bonded conjugates, to give a definitive baseline for comparison to the activity of targeted toxins from the BETA product pipeline as well as others. BIgG-SAP serves as a control for bonded mouse conjugates utilizing biotin-streptavidin technology.

Specificity and Preparation:

This control conjugate has no known specificity. BIgG-SAP mouse is a bonded conjugate of biotinylated mouse IgG and the secondary conjugate Streptavidin-ZAP (IT-27) containing and the ribosome-inactivating protein, saporin.

Usage and Storage:

BIgG-SAP serves as a control for bonded mouse conjugates utilizing biotin-streptavidin technology. **There may be lot-to-lot variation in material; working dilutions must be determined by end user. If this is a new lot, you must assess the proper working dilution before beginning a full experimental protocol.** Gently spin down material before use; 5-10 seconds in a microfuge should be adequate. The material should be stored at -20°C in undiluted aliquots. Material should be aliquoted to a convenient volume and quantity to avoid repeated freezing and thawing that can damage the protein content. Under these conditions, the material has a very stable shelf-life. Thawing should be done at room temperature or on ice. The thawed solution should remain on ice until use. The material can be handled safely using normal laboratory precautions. Do not use a reducing agent (such as dithiothreitol, beta-mercaptoethanol or ascorbic acid) with this material. It will inactivate the toxin.
For disposal: autoclave, or expose to 0.2 M NaOH, materials that come into contact with the toxin.

Safety:

Good laboratory technique must be employed for safe handling of this product.

This requires observation of the following practices:

1. Wear appropriate laboratory attire, including lab coat, gloves and safety glasses.
2. Do not pipet by mouth, inhale, ingest or allow product to come into contact with open wounds. Wash thoroughly any part of the body which comes into contact with the product.
3. Avoid accidental autoinjection by exercising extreme care when handling in conjunction with any injection device.
4. This product is intended for research use by qualified personnel only. It is not intended for use in humans or as a diagnostic agent. Advanced Targeting Systems is not liable for any damages resulting from the misuse or handling of this product.

To view protocol(s) for this and other products please visit: www.ATSBio.com/support/protocols



6352 CORTE DEL ABETO, STE B
CARLSBAD, CA 92011 USA
01.858.642.1988 • WWW.ATSBIO.COM

BIgG-SAP Mouse
NON-TARGETED SAPORIN CONTROL MOLECULE