

**Antibody to Mouse CD90 Thy-1**
RAT MONOCLONAL

Catalog Number: PRP-142ANT
Quantity: 500 micrograms
Format: Lyophilized
Host: Rat
Isotype: IgG1
Clone: mThy-1
Immunogen: Purified mouse LN T cells

Background: CD-90 plays a role in cell-cell or cell-ligand interactions during synaptogenesis and other events in the brain. Thy-1 is the designation for a major cell surface glycoprotein characteristic to T cells. The Thy-1 glycoproteins are constituents of thymocytes and neurons and probably are involved in cell-cell interactions. Human Thy-1 homolog maps to chromosome 11 because that is where T3D is found, in the mouse T3D and Thy-1 map to chromosome 9 along with certain other loci that are on human 11q.

Specificity & Preparation: Rat anti-mouse CD90 Thy-1 is purified by ion exchange column. Protein concentration is 1 mg/ml in PBS (after reconstitution).

Usage: Reported to be effective as a cytotoxic and staining antibody. For staining, use 10 μ l/10⁶ cells. This antibody is extremely efficient for killing *in vivo*. Exact titer for cytotoxicity should be determined by the investigator.

Storage: Reconstitute with H₂O. Gently spin down material 5-10 seconds in a microfuge before use. Store lyophilized material at 4°C. After reconstitution, if not intended for use within a month, aliquot and store at -20°C. Material is stable for two years lyophilized, one month in solution at 4°C.

To view protocol(s) for this and other products please visit: www.ATSBio.com/library/protocols