

**Biotin-labeled Antibody to Human CD4**
MOUSE MONOCLONAL

Catalog Number: PRP-167ANT
Quantity: 500 micrograms
Format: Lyophilized
Host: Mouse
Isotype: IgG2b
Clone: hCD4
Immunogen: Purified human PBL CD4+ T cells

Background: CD4 is a cell-surface glycoprotein found on mature helper T-cells, immature thymocytes, monocytes, and macrophages. Some cytotoxic T-cells express CD4 as well. The T-cell population in normal blood is approximately 65% CD4+. A mature T-cell expresses either CD4 or CD8, but not both. During an early stage of development T-cells have both CD4 and CD8 receptors but are differentiated in the thymus.

Specificity & Preparation: Mouse anti-human CD4 is conjugated to biotin and purified by ion exchange chromatography. The protein concentration is 1 mg/ml in PBS after reconstitution.

Usage: Reported to be effective as a blocking and staining antibody. Blocks binding of HIV to CD4. For staining, use 10 μ l/10⁶ cells. The titer for blocking T-cell activation should be determined by the investigator.

Storage: Reconstitute with H₂O. Gently spin down material 5-10 seconds in a microfuge before use. Store lyophilized material at 4°C. After reconstitution, if not intended for use within a month, aliquot and store at -20°C. Material is stable for two years lyophilized, one month in solution at 4°C.

To view protocol(s) for this and other products please visit: www.ATSBio.com/library/protocols