



6352 CORTE DEL ABETO, STE B
CARLSBAD, CA 92011 USA
01.858.642.1988 • WWW.ATSBIO.COM

RAB5A Human
RECOMBINANT PROTEIN

Catalog Number: PRP-590PRO
Quantity: 5 micrograms, 20 micrograms, 1 milligram
Format: Sterile-filtered colorless solution containing 20mM Tris, pH-8
Host: *E. coli*

Background:

RAB5A is a member of RAS oncogene family. The small GTPase, RAB5A is a crucial regulatory component of the endocytic pathway. Activation of RAB5A is mediated by GDP-GTP exchange factors (GEFs) that generate the Rab5-GTP complex). RAB5A is required for the fusion of plasma membranes and early endosomes.

Specificity and Preparation:

RAB5A Human Recombinant produced in *E.Coli* is a single, non-glycosylated polypeptide chain containing 215 amino acids and having a molecular mass of 23kDa. Greater than 95.0% purity as determined by analysis by RP-HPLC and SDS-PAGE.

Amino Acid Sequence:

MASRGATRPN GPNTGNKICQ FKLVLGESA VGKSSLVLRV VKGQFHEFQE
STIGAAFLTQTVCLDDTTVK FEIWDTAGQE RYHSLAPMYR RGAQAIVVY DITNEESFAR
AKNWKELQRQASPNIVIAL SGNKADLANK RAVDFQEAQS YADDNSLLFM ETSKTSMTNV
NEIFMAIAKKLPKNPQNPQ ANSARGRGVD LTPPTQPTRN QCCSN.

Usage and Storage:

RAB5A although stable at 4°C for 4 weeks, should be stored desiccated below -18°C. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid repeated freezing and thawing. Gently spin down material before use; 5-10 seconds in a microfuge should be adequate.

To view protocol(s) for this and other products please visit: www.ATSBio.com/support/protocols